

Image-Pro v11

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	Pricing (USD)
All Products have the option of being delivered with a USB License Dongle, listed as Physical (P) or delivered Electronically via E-Delivery (E) to be linked to the Computer for the same price .				

BUNDLES

Bundles - General				
100-107	Image-Pro - Capture Bundle	Capture images & videos from 150+ devices and then process, adjust, and view images with this powerful and streamlined capture bundle. INCLUDES: Base Product (100-001), 2D Capture Module (300-016)	P -or- E	\$900
100-101	Image-Pro - Basic Bundle	Extract, compile, review, and report on precise manual measurement data with entry-level Basic bundle. INCLUDES: Base Product (100-001), 2D Measurements Module (300-011)	P -or- E	\$1,530
100-102	Image-Pro - 2D Analysis Bundle	Extract, segment, classify, and analyze precise manual and automated measurement data with the 2D Analysis bundle. NOTE: This bundle is equivalent to and replaces the "Image-Pro Single User License (USB)" (100-002) INCLUDES: Base Product (100-001), 2D Measurements Module (300-011), 2D Automated Analysis Module (300-012)	P -or- E	\$5,750
100-103	Image-Pro - 2D/3D Analysis Bundle	Extract, segment, classify, and analyze precise manual and automated measurement data from 2D or 3D images with the 2D/3D Analysis bundle. INCLUDES: Base Product (100-001), 2D Measurements Module (300-011), 2D Automated Analysis Module (300-012), 3D Visualization Module (300-019), & 3D Measure & Analyze Extension (300-020)	P -or- E	\$12,500
Bundles - Life Science				
100-108	Image-Pro - Multi-Channel Capture Bundle	Expand any microscope into a multi-channel fluorescence capture system with this bundle. Control lightsource devices, capture mixed sets of fluorescence and brightfield, and deconvolve using AutoQuant Real-Time Deconvolution, resulting in a crisp, high fidelity fluorescent image set. INCLUDES: Base Product (100-001), 2D Capture - Module (300-016), Multi-Channel Capture - Extension (300-018), & Real-Time Deconvolution - Extension (300-017)	P -or- E	\$1,620
100-104	Image-Pro - Deconvolution Station Bundle	Import, deconvolve, visualize, and generate animations in a standalone bundle designed to process and visualize fluorescent 3D image stacks. INCLUDES: Base Product (100-001), 3D Visualization Module (300-019), & AutoQuant Deconvolution Module (300-022)	P -or- E	\$4,770
100-105	Image-Pro - 2D Cell Analysis Bundle	Extract, segment, classify, analyze, and report precise 2D manual and automated cell measurement data with the 2D Cell Analysis bundle. INCLUDES: Base Product (100-001), 2D Measurements Module (300-011), 2D Automated Analysis Module (300-012), Cell Biology Protocol Extension (300-013), Cell Biology Plus Protocol Extensions (300-014), & AutoQuant Deconvolution Module (300-022)	P -or- E	\$9,440
Bundles - Industrial				
100-106	Image-Pro - 2D Materials Analysis Bundle	Extract, segment, classify, analyze, and report precise 2D manual and automated materials measurement data with the 2D Materials Analysis bundle. INCLUDES: Base Product (100-001), 2D Measurements Module (300-011), 2D Automated Analysis Module (300-012), Materials Protocol Extension (300-015)	P -or- E	\$6,336

Image-Pro v11

INDIVIDUAL

Base Product				
100-001	Image-Pro - Base	The Image-Pro Base is required for all modules and extensions. Includes File Management, 2D Viewer, Process/Adjust, Usability, Data Management, Reports, Auditing, Scripting LINK TARGETS: [Windows PC, Green USB dongle, or License Server Module (sold separately); Purple USB dongles not supported]	P -or- E	\$300
Modules & Extensions				
300-011	2D Measurements - Module	Extract precise 2D manual measurement data using semi-automated tools for fast and accurate results. Includes Direct & Assisted Measurements, Auto Shape & Edge Detection, Classification, & Color Measurement. REQUIRES: [100-001 - Image-Pro - Base]	P -or- E	\$1,400
300-012	2D Automated Analysis - Module	Tackle complex analysis across multiple images/channels with Protocols or Free-Form analysis with ML-based tools. Includes Smart Segmentation, Classification, 2D Tracking, & the Essentials Protocol Collection with 3 Protocols. REQUIRES: [100-001 - Image-Pro - Base]	P -or- E	\$4,690
300-013	Protocols - Cell Biology - Extension	Analyze complex Cell Biology-based measurements with an efficient and repeatable process. Includes 10 Protocols: Apoptosis, Autophagy, Cell Count, Cell Morphology, Cell Proliferation, Confluence, Lipid Droplet, Live/Dead, Ring, & Transfection. REQUIRES: [300-012 - 2D Automated Analysis - Module]	P -or- E	\$500
300-014	Protocols - Cell Biology Plus - Extension	Analyze complex Cell Biology-based measurements with an efficient and repeatable process. Includes 5 Protocols: Angiogenesis, Colocalization, Neurite Outgrowth, Translocation, & Wound Healing. REQUIRES: [300-012 - 2D Automated Analysis - Module]	P -or- E	\$600
300-015	Protocols - Materials - Extension	Analyze complex Materials-based measurements with an efficient and repeatable process. Includes 10 Protocols: Composite Materials, Layers Concentric, Layers Curved, Layers Parallel, Fiber Separation, Fiber Thickness, Particles, Particles Phase, Phase, & Pores. REQUIRES: [300-012 - 2D Automated Analysis - Module]	P -or- E	\$650
300-016	2D Capture - Module	Capture images & videos from 150+ devices, perform live measurement, and live processing with ease. Includes the Live Image Workspace, Instant & HDR Capture, Recording of Time-Lapse Events, and basic Presets, Enhancements, and Preview/Acquire settings. REQUIRES: [100-001 - Image-Pro - Base]	P -or- E	\$700
300-017	Real-Time Deconvolution - Extension	Elevate any capture system with Real-Time Deconvolution to deliver crisp, high fidelity, fluorescent images during preview and capture. Includes 2D iterative Deconvolution & Wiener Filter algorithms. REQUIRES: [300-016 - 2D Capture - Module]	P -or- E	\$400
300-018	Multi-Channel Capture - Extension	Expand any microscope with Multi-Channel Capture making it capable of capturing multiple channels into an image set. Includes LightSource Device Control, Acquisition of mixed sets with Fluorescence and Brightfield, and Auto & Manual capture modes. REQUIRES: [300-016 - 2D Capture - Module]	P -or- E	\$400

Image-Pro v11

300-019	3D Visualization - Module	<p>Visualize 3D image stacks with full control over appearance for the generating of stunning animations with complex movements & renderings. Includes 3D viewer with advanced controls, Min/Max projections, clipping planes, slicer viewing, & Movie Maker.</p> <p>REQUIRES: [100-001 - Image-Pro - Base]</p>	P -or- E	\$2,000
300-020	3D Measure & Analyze - Extension	<p>Measure & Analyze with familiar tools across multiple dimensions. Includes iso-surface based object segmentation, 3d volume measurements, guided segmentation, 3D colocalization, and manual point-to-point distances, angles, and more in 3D space.</p> <p>REQUIRES: [300-019 - 3D Visualization - Module]</p>	P -or- E	\$5,500
300-021	3D Big Data Visualization - Extension	<p>Convert large datasets captured by high-volume imaging techniques to multi-resolution pyramidal formats for quick and easy visualization. Includes background file conversion, terabyte-sized file support, and multi-resolution file visualization.</p> <p>REQUIRES: [300-019 - 3D Visualization - Module]</p>	P -or- E	\$2,000
300-022	AutoQuant Deconvolution - Module	<p>Restore image fidelity, enhance quality, and reverse distortions introduced during the image acquisition process. Includes 2D & 3D Algorithms, Theoretical & Measured PSF models, Spherical Abberation correction, GPU acceleration, and Batch processing.</p> <p>REQUIRES: [100-001 - Image-Pro - Base]</p>	P -or- E	\$3,000
300-004	License Server - Module	<p>Distribute software licenses to teams with a stand-alone License Server and remove USB dongle and Node-Locked licenses from the network. Includes floating or reserved to workstation options, flexible bundle assignment, auto-license return after user time outs, options for VPN and Virtual Machine support.</p> <p>REQUIRES: [100-001 - Image-Pro - Base]</p> <p>LINK TARGETS: [Windows, Windows Server, Linux, VM's ; See IT Checklist for complete list]</p>	P -or- E	\$2,000

Image-Pro v11

Upgrades				
100-016	Image-Pro - Upgrade from Basic & Insight Products	<p>Upgrade to the latest version of Image-Pro from any previous version of a "Basic" or "Insight" branded product to receive a product with similar features and configuration. If the previous configuration included a Capture Module, the upgrade to the latest version will also include the new equivalent module at no additional charge.</p> <p>REQUIRES: Valid Serial Number</p> <p>INCLUDES: Image-Pro - Basic Bundle (100-101), Optional: 2D Capture - Module (300-016)</p>	P -or- E	\$400
100-020	Image-Pro - Upgrade from Plus & Analyzer v6.3 and older Products	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Plus" & "Analyzer" v6.3 and older to receive a product with similar features and configuration. If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent module at no additional charge.</p> <p>REQUIRES: Valid Serial Number</p> <p>INCLUDES: Image-Pro - 2D Analysis Bundle (100-102), Optional: 2D Capture - Module (300-016), Optional: License Server - Module (300-004)</p>	P -or- E	\$5,290
100-017	Image-Pro - Upgrade from Plus & Analyzer v7 Products	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Plus" or "Analyzer" v7 to receive a product with similar features and configuration. If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent module at no additional charge.</p> <p>REQUIRES: Valid Serial Number</p> <p>INCLUDES: Image-Pro - 2D Analysis Bundle (100-102), Optional: 2D Capture - Module (300-016), Optional: License Server - Module (300-004)</p>	P -or- E	\$4,200
100-018	Image-Pro - Upgrade from v9, 10, and AutoQuant Products	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Premier" v9 or Image-Pro v10 or even any version of AutoQuant to receive a product with similar features and configuration. If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent module at no additional charge.</p> <p>REQUIRES: Valid Serial Number</p> <p>INCLUDES: Image-Pro - Deconvolution Station Bundle -OR- Image-Pro - 2D Analysis Bundle (100-102), Optional: 2D Capture - Module (300-016), Optional: License Server - Module (300-004)</p>	P -or- E	\$2,350
100-019	Image-Pro - Upgrade from v9 or v10 Products - with 3D module	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Premier 3D" v9 or Image-Pro v10 with the "3D module" to receive a product with similar features and configuration. If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent module at no additional charge.</p> <p>REQUIRES: Valid Serial Number</p> <p>INCLUDES: Image-Pro - 2D/3D Analysis Bundle (100-103), Optional: 3D Big Data Visualization - Extension (300-021), License Server - Module (300-004)</p>	P -or- E	\$3,350

Image-Pro v11

Support & Accessories				
700-002	Upgrade Protection Plan	Upgrade Protection ensures your product is always up-to-date with the latest features. While under coverage, all upgrades to features within your product configuration will be provided at no additional cost. Annual Renewal required. REQUIRES: [100-001 - Image-Pro - Base]	E	\$945
700-016	USB Dongle	A single secure USB dongle (green) with license linking capability for Image-Pro. Use the License Manager to select this USB dongle as a target for linking an Image-Pro v11 or higher license.	P	\$700
700-017	USB Dongle - Demo NFR - Black	A single secure USB dongle (black) that can only be linked to a Not For Resale (NFR) demo license. Use this dongle for Demo Purposes ONLY. Use the License Manager to select this USB dongle as a target for linking an Image-Pro v11 or higher NFR Demo License.	P	\$150

Image-Pro Plus & Analyzer

Product Code	Product Name	Description	Physical / E-Delivery	Pricing (USD)
NEW				
Single Permanent Licenses				
102-001	Image Pro Plus v7.1 - Single User License (EN, CN, JP) - USB Dongle	Installation USB Flash Drive & USB License Dongle	Physical	\$7,268
102-005	Image-Pro Analyzer v7.1 - Single User License (EN) - USB Dongle	Installation USB Flash Drive & USB License Dongle	Physical	\$6,156
Upgrades				
102-008	Image-Pro Plus v7.1 (EN, CN, JP) - Upgrade from v6.3 thru v5.0	An Upgrade Authorization Code (UAC) will be provided along with a link to download the installer. UAC can be used to upgrade an existing USB license dongle. If no USB dongle is available, use 700-015 to purchase new USB dongle. REQUIRES: Valid Serial Number	E-Delivery	\$5,390
102-009	Image-Pro Analyzer v7.1 (EN) - Upgrade from v6.3 thru v5.0	An Upgrade Authorization Code (UAC) will be provided along with a link to download the installer. UAC can be used to upgrade an existing USB license dongle. If no USB dongle is available, use 700-015 to purchase new USB dongle. REQUIRES: Valid Serial Number	E-Delivery	\$4,755
102-010	Image-Pro Plus v7.1 (EN, CN, JP) - Upgrade from v7.0	An Upgrade Authorization Code (UAC) will be provided along with a link to download the installer. UAC can be used to upgrade an existing USB license dongle. If no USB dongle is available, use 700-015 to purchase new USB dongle. REQUIRES: Valid Serial Number	E-Delivery	\$1,500
102-011	Image-Pro Analyzer v7.1 (EN) - Upgrade from v7.0	An Upgrade Authorization Code (UAC) will be provided along with a link to download the installer. UAC can be used to upgrade an existing USB license dongle. If no USB dongle is available, use 700-015 to purchase new USB dongle. REQUIRES: Valid Serial Number	E-Delivery	\$1,275
Accessories				
700-015	Replacement USB Dongle for Image-Pro v7 - v10, & AutoQuant X3	A single secure USB dongle	Physical	\$675