




Image-Pro v11.1 – NFR License

Registration, Downloading & Linking

STEP #1: Download & Install Image-Pro v11.1 & AI Neural Engine



<https://my.mediacy.com/support/updates/image-pro-11-1-0-426>



STEP #2: Open Image-Pro & Login Using License Manager

Note: Use same login credentials as my.MediaCy.com

If you did not receive an account registration email or need an account invitation, please contact Media Cybernetics

License Manager

Internet Connection: Online Account Status: Sign In

License Manager

Sign In
with myMediaCy Account

• Full access to licenses in your account

Username
tmoyer@mediacy.com

Password Show

Sign in

OR

Continue
with Unlocking Codes

This Computer
Installed Software
My Unlinked
My Linked L
Software

Target
Link
Action
Unlink



STEP #3: Select an Available Image-Pro – NFR License

Note:

Choose an available license under the “My Unlinked Licenses” menu

The screenshot shows the License Manager interface with the following sections:

- This Computer:** A table showing installed software. One entry is 'Image-Pro' with version 11.1 and status 'Linked'.
- My Unlinked Licenses:** A table with columns: Software Licenses, ULC, Version, Status, and Action. It lists three 'Image-Pro - NFR Demo (195 days left)' licenses with ULCs ...6K6TQ, ...4KVCT, and ...WC9R9. The status for the first two is 'Unlinked', and the third is also 'Unlinked' (highlighted with a blue box). There is a 'Link' button at the bottom right of this section.
- My Linked Licenses:** A table with columns: Software Licenses, ULC, Version, Status, and Action. It lists three 'Image-Pro - NFR Demo (195 days left)' licenses with ULCs ...9FYCP, ...3XYH6, and ...WC9R9. The status for the first is 'Linked', the second is 'Update', and the third is 'Linked'. There is an 'Unlink' button at the bottom right of this section.

Note:

Make sure the license status is “unlinked”



STEP #4: Click “Link”

The screenshot shows the License Manager application window. At the top, it displays 'Internet Connection: Online' and 'Account Status: Tom Moyer'. Below this are navigation tabs for 'My Demo Products', 'License Server', and 'My Users'. The main content area is divided into three sections: 'This Computer', 'My Unlinked Licenses', and 'My Linked Licenses'.

This Computer

Installed Software	Version	Status	Target
Image-Pro	11.1	✓ Linked	

My Unlinked Licenses

<input checked="" type="checkbox"/>	Software Licenses	ULC	Version	Status	Action
<input checked="" type="checkbox"/>	Image-Pro - NFR Demo (195 days left)	...9FYCP	11.1	• Unlinked	Link
<input type="checkbox"/>	Image-Pro - NFR Demo (195 days left)	...6K6TQ	11.1	• Unlinked	
<input type="checkbox"/>	Image-Pro - NFR Demo (195 days left)	...4KVCT	11.1	• Unlinked	

My Linked Licenses

<input type="checkbox"/>	Software Licenses	ULC	Version	Status	Action
<input type="checkbox"/>	Image-Pro - NFR Demo (195 days left)	...3XYH6	11.1	Update	
<input type="checkbox"/>	Image-Pro - NFR Demo (195 days left)	...WC9R9	11.1	Linked	
<input type="checkbox"/>	Image-Pro - NFR				Unlink

A red arrow points to the 'Link' button in the 'My Unlinked Licenses' table.



STEP #5: Choose USB or PC Linking Option

The screenshot shows the License Manager application window. At the top, it displays 'Internet Connection: Online' and 'Account Status: Tom Moyer'. Below this are navigation tabs for 'My Demo Products', 'License Server', and 'My Users'. The main content area is divided into three sections:

- This Computer:** A table showing installed software. One entry is 'Image-Pro' with version 11.1 and status 'Linked'.
- My Unlinked Licenses:** A table of unlinked licenses. A 'Choose Link Target' dialog box is open over this section. The dialog has two radio buttons: 'This Computer' (selected) and 'USB Dongle'. The 'USB Dongle' option has a note: 'USB Dongle is missing. Plug in to detect.' Below the dialog, a table shows unlinked licenses, including 'Image-Pro - NFR Demo (195 days left)' with ULC '...4KVCT' and version 11.1. A 'Link' button is visible at the bottom right of this section.
- My Linked Licenses:** A table of linked licenses. It shows 'Image-Pro - NFR Demo (195 days left)' with ULC '...3XYH6' and version 11.1, and another entry with ULC '...WC9R9' and version 11.1. An 'Unlink' button is at the bottom right.



STEP #6: Close License Manager & Run Image-Pro v11.1



Please contact Media Cybernetics with any questions

Tom Moyer
Media Cybernetics | Sales Manager

1700 Rockville Pike, Suite 240,
Rockville, MD 20852

Phone: +1-301-495-3305 Ext. 223
tmoyer@mediacy.com

