

Image-Pro v11.1

Table of Contents

PACKAGES

PACKAGES - OVERVIEW

[Worksheet - Packages Overview](#)

IMAGE PRO - AI DEEP LEARNING PACKAGES

LIFE SCIENCE PACKAGES | MATERIALS PACKAGES | GENERAL PACKAGES

[Worksheet - AI Packages](#)

IMAGE PRO - MACHINE LEARNING PACKAGES

LIFE SCIENCE PACKAGES | MATERIALS PACKAGES | GENERAL PACKAGES

[Worksheet - ML Packages](#)

IMAGE PRO - AUTOQUANT PACKAGES

AUTOQUANT PACKAGES

[Worksheet - AutoQuant Packages](#)

IMAGE PRO - CAPTURE PACKAGES

LIFE SCIENCE CAPTURE PACKAGES | GENERAL CAPTURE PACKAGES

[Worksheet - Capture Packages](#)

SUCCESS PLANS

SUCCESS PLANS - OVERVIEW

[Worksheet - Success Plans Overview](#)

STANDARD SUCCESS PLANS | PROFESSIONAL SUCCESS PLANS | ULTIMATE SUCCESS PLANS

[Worksheet - Success Plans](#)

UPGRADES

UPGRADES FROM LEGACY VERSIONS | UPGRADES FROM v9, v10 & AUTOQUANT | IMAGE PRO - UPGRADES FROM v11.0

[Worksheet - Upgrades](#)

COMPLETE PRICE LIST

PACKAGES | INDIVIDUAL PRODUCTS | UPGRADES | ACCESSORIES | SERVICE PLANS

[Worksheet - Complete Price List](#)

Image-Pro v11.1

PACKAGES - OVERVIEW

Name		General					Life Science					Materials			
		Basic	Capture	2D Analysis	3D Analysis	2D/3D Analysis	AI Complete Analysis (2D/3D)	Multi-Channel Capture	Deconvolution Station	ML Life Science (2D)	AI Life Science (2D)	AI Life Science (2D/3D)	ML Mat Analysis (2D)	AI Materials (2D)	AI Materials (2D/3D)
Part #s		100-205	100-204	100-206	100-219	100-207	100-220	100-208	100-209	100-210	100-211	100-221	100-212	100-213	100-222
Package List Price (USD)		\$1,575	\$945	\$5,795	\$7,065	\$12,545	\$19,520	\$1,665	\$4,905	\$6,875	\$11,060	\$17,810	\$6,425	\$10,655	\$17,405
Image-Pro	100-201	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
2D Measurements	300-024	✓		✓		✓	✓			✓	✓	✓	✓	✓	✓
2D Automated Analysis	300-025			✓		✓	✓			✓	✓	✓	✓	✓	✓
Protocols - Cell Biology	300-026						✓			✓	✓	✓			
Protocols - Cell Biology Plus	300-027						✓			✓	✓	✓			
Protocols - Cell Biology EM	300-032						✓			✓	✓				
Protocols - Materials	300-028						✓					✓	✓	✓	
Protocols - Materials EM	300-033						✓						✓	✓	
AI Deep Learning	300-029						✓				✓	✓		✓	✓
Pre-Trained Models - Life Science	300-030						✓				✓	✓			
Pre-Trained Models - Materials	300-031						✓						✓	✓	
2D Capture	300-034		✓					✓							
Real-Time Deconvolution	300-035							✓							
Multi-Channel Capture	300-036							✓							
3D Visualization	300-037				✓	✓	✓		✓			✓			✓
3D Measure & Analyze	300-038				✓	✓	✓					✓			✓
3D Big Data Visualization	300-039														
AutoQuant Deconvolution	300-040								✓						
License Server	300-023														

Image-Pro v11.1

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
--------------	--------------	-------------	----------------------------------	---------------------

All Products have the option of being delivered with a USB License Dongle, listed as **Physical (P)** or delivered Electronically via **E-Delivery (E)** to be linked to the Computer for the **same price**.

IMAGE PRO - AI DEEP LEARNING PACKAGES

IMAGE PRO - AI LIFE SCIENCE PACKAGES

100-211	Image-Pro v11.1 - AI Life Science Package (2D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Life Science applications with the AI Life Science Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Protocols - Cell Biology, Protocols - Cell Biology Plus, Protocols - Cell Biology EM, AI Deep Learning Extension, Pre-Trained Models - Life Science	P -or- E	\$11,060
100-221	Image-Pro v11.1 - AI Life Science Package (2D/3D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Life Science applications and is combined with 3D Visualization and Analysis in the AI Life Science Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Protocols - Cell Biology, Protocols - Cell Biology Plus, Protocols - Cell Biology EM, AI Deep Learning Extension, Pre-Trained Models - Life Science, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$17,810

IMAGE PRO - AI MATERIALS PACKAGES

100-213	Image-Pro v11.1 - AI Materials Package (2D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Materials applications with the AI Materials Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlysis Module, Materials Protocols, & Materials EM Protocol Extension, AI Deep Learning Module, Pre-Trained Materials Models Extension	P -or- E	\$10,655
100-222	Image-Pro v11.1 - AI Materials Package (2D/3D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Materials applications and is combined with 3D Visualization and Analysis in the AI Materials Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlysis Module, Materials Protocols, & Materials EM Protocols Extension, AI Deep Learning Module, Pre-Trained Materials Models Extension, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$17,405

IMAGE PRO - AI GENERAL PACKAGES

100-220	Image-Pro v11.1 - AI Complete Analysis Package (2D/3D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D Life Science and Materials applications and is combined with 3D Visualization and Analysis in the AI Complete Analysis Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Protocols - Cell Biology, Protocols - Cell Biology Plus, Protocols - Cell Biology EM, Protocols - Materials, Proocols - Materials EM, AI Deep Learning Extension, Pre-Trained Models - Life Science, Pre-Trained Models - Materials, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$19,520
---------	--	--	----------	----------

Image-Pro v11.1

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
--------------	--------------	-------------	-------------------------------	------------------

All Products have the option of being delivered with a USB License Dongle, listed as **Physical (P)** or delivered Electronically via **E-Delivery (E)** to be linked to the Computer for the **same price**.

IMAGE PRO - MACHINE LEARNING PACKAGES

IMAGE PRO - ML LIFE SCIENCE PACKAGES				
100-210	Image-Pro v11.1 - Life Science Package (2D)	Machine Learning powered segmentation and classification enables precise analysis and automated measurement from 2D images with the ML Life Science Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Cell Biology, Cell Biology Plus Protocol Extensions	P -or- E	\$6,875

IMAGE PRO - ML MATERIALS PACKAGES				
100-212	Image-Pro v11.1 - Materials Package (2D)	Machine Learning powered segmentation and classification enables precise analysis and automated measurement from 2D images with the ML Materials Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlysis Module, Materials Protocol Extension	P -or- E	\$6,425

IMAGE PRO - ML GENERAL PACKAGES				
100-205	Image-Pro v11.1 - Basic Package	Extract, compile, review, and report on precise manual measurement data with entry-level Basic bundle. INCLUDES: Base Product, 2D Measurements Module	P -or- E	\$1,575
100-206	Image-Pro v11.1 - 2D Analysis Package	Machine Learning powered segmentation and classification enables precise analysis and automated measurement with the 2D Analysis package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module	P -or- E	\$5,795
100-219	Image-Pro v11.1 - 3D Analysis Package	Visualize stunning 3D views with robust 3D controls and perform precise manual 3D measurements and automated 3D segmentation with the 3D Analysis package. INCLUDES: Base Product, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$7,065
100-207	Image-Pro v11.1 - 2D/3D Analysis Package	Machine Learning powered segmentation and classification enables precise analysis and automated measurement from 2D or 3D images with the 2D/3D Analysis package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$12,545

Image-Pro v11.1

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
--------------	--------------	-------------	----------------------------------	---------------------

All Products have the option of being delivered with a USB License Dongle, listed as **Physical (P)** or delivered Electronically via **E-Delivery (E)** to be linked to the Computer for the **same price**.

IMAGE PRO - DECONVOLUTION PACKAGES

IMAGE PRO - AUTOQUANT PACKAGES				
100-209	Image-Pro v11.1 - Deconvolution Station Package	Import, deconvolve, visualize, and generate animations in a standalone bundle designed to process and visualize fluorescent 3D image stacks. INCLUDES: Base Product, 3D Visualization Module, & AutoQuant Deconvolution Module	P -or- E	\$4,905

Image-Pro v11.1

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
--------------	--------------	-------------	----------------------------------	---------------------

All Products have the option of being delivered with a USB License Dongle, listed as **Physical (P)** or delivered Electronically via **E-Delivery (E)** to be linked to the Computer for the **same price**.

IMAGE PRO - CAPTURE PACKAGES

IMAGE PRO - LIFE SCIENCE CAPTURE PACKAGES				
100-208	Image-Pro v11.1 - Multi-Channel Capture Package	Expand any microscope into a multi-channel fluorescence capture system with this package. Control lightsource devices, capture mixed sets of fluorescence and brightfield, and deconvolve using AutoQuant Real-Time Deconvolution, resulting in a crisp, high fidelity fluorescent image set. INCLUDES: Base Product, 2D Capture - Module, Multi-Channel Capture - Extension, & Real-Time Deconvolution - Extension	P -or- E	\$1,665

IMAGE PRO - GENERAL CAPTURE PACKAGES				
100-204	Image-Pro v11.1 - Capture Module	Capture images & videos from 150+ devices and then process, adjust, and view images with this powerful and streamlined capture package. INCLUDES: Base Product, 2D Capture Module	P -or- E	\$945

Image-Pro v11.1

SUCCESS PLANS - OVERVIEW

SUCCESS PLANS - 1 YEAR			
	STANDARD	PROFESSIONAL	ULTIMATE
Success Plan - Single License <i>Needs to cover all licenses at facility</i>	700-018	700-020	700-022
	\$1,450	\$1,850	\$3,895
Success Plan - Additional Licenses <i>Add Per Additional License</i>	700-019	700-021	700-023
	\$150	\$1,095	\$1,250

SUCCESS PLANS - 3 YEARS (20% DISCOUNT)			
	STANDARD	PROFESSIONAL	ULTIMATE
Success Plan - Single License <i>Needs to cover all licenses at facility</i>	700-024	700-026	700-028
	\$3,480	\$4,440	\$9,348
Success Plan - Additional Licenses <i>Add Per Additional License</i>	700-025	700-027	700-029
	\$360	\$2,628	\$3,000

Image-Pro v11.1

SUCCESS PLANS - FEATURES

	STANDARD	PROFESSIONAL	ULTIMATE
Software Onboarding	✓	✓	✓
Technical Support	✓	✓	✓
Software Upgrade Protection		✓	✓
New Protocol Pass		✓	✓
New AI Model Pass		✓	✓
Uptime Guarentee			✓
License Server Support			✓
Software Training Hours			✓
Application Service Hours			✓
AI Model Training-as-a-Service (option)			✓

Image-Pro v11.1

SUCCESS PLANS - FEATURE DESCRIPTIONS

Software Onboarding	Onboarding includes a single session to complete software installation, configuration, and basic usage orientation. Services are provided by Media Cybernetics professional support staff by remote meeting to help ensure your business can confidently navigate and leverage the software's features, driving efficiency and immediate value.
Technical Support	Technical Support includes professional and timely assistance to troubleshoot issues, help with proper software configuration, set up 3rd party hardware (i.e. cameras), and answer usage questions to help you get the most out of your investment.
Software Upgrade Protection	Upgrade Protection ensures your product is always up-to-date with the latest features. While under coverage, all new version upgrades will be provided at no additional cost within the existing product configuration. To receive modules outside of the existing product configuration, a separate purchase is required.
New Protocol Pass	The New Protocol Pass grants continuous access to the latest image analysis protocols for your Image-Pro software and licensed Protocol Collections during the term of your plan. Pass holders receive new protocols as they are developed, enhancing software capabilities with cutting-edge tools and techniques.
New AI Model Pass	The New AI Model Pass provides continuous access to new pre-trained AI models for your Image-Pro software and licensed AI Model Extensions during the term of your plan. Pass holders stay updated with the latest AI model advancements, maximizing the potential of your software's functionality.
Uptime Guarantee	The Uptime Guarantee provides a minimum of three 30-day Flex licenses of Image-Pro v11.0 and newer. Additional Flex licenses are available upon request. All Flex licenses will expire at the end of their term, regardless of usage. These licenses can be used when permanent licenses are unavailable due to computer maintenance or to accommodate fluctuating software license needs throughout the year.
License Server Support	A Technical Support Single Incident (for License Server issues) gives access to the Media Cybernetics professional support staff to tackle a single issue related to the License Server from start to finish. The issue can include troubleshooting a License Server problem, helping with proper License Server configuration, working with IT to re-deploy a server, or simply answering a usage question.
Software Training Hours	Up to 2 hours of online training delivered by specialists who know your software the best. This training offers an overview of the Image-Pro product with a tailored agenda to address your specific needs.
Application Service Hours	Up to 2 hours of expert guidance on deploying Image-Pro software. Advice aims to provide best practices, custom solutions as needed, and process optimizations to streamline your workflow and maximize your investment in Image-Pro. Does not include software development or macro writing. *Hours can be used for AI Model Training-as-a-Service.
AI Model Training-as-a-Service (option)	Expert training of AI Deep Learning CNN models for microscopic image analysis, including consultation, image annotation, and model training. The service ensures high-quality results, saving time and resources, and includes a detailed report and marked-up images. *Can be used in place of Application Service Hours.

Image-Pro v11.1

Product Code	Product Name	Description	Product Type	List Price (USD)
--------------	--------------	-------------	--------------	------------------

SUCCESS PLANS

STANDARD SUCCESS PLANS				
700-018	Success Plan - Standard (Single License) - 1 Year	<p>This plan provides a single onboarding session and continuous software support for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-019 - Standard - Additional License]</p>	Subscription	\$1,450
700-019	Success Plan - Standard (Additional License) - 1 Year	<p>Add an additional license to your Plan and receive continuous software support for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>REQUIRES: [700-018 - Standard (Single License)]</p>	Subscription	\$150
700-024	Success Plan - Standard (Single License) - 3 Years	<p>This plan provides a single onboarding session and continuous software support for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-025 - Standard - Additional License]</p>	Subscription	\$3,480
700-025	Success Plan - Standard (Additional License) - 3 Years	<p>Add an additional license to your Plan and receive continuous software support for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>REQUIRES: [700-024 - Standard (Single License)]</p>	Subscription	\$360

Image-Pro v11.1

PROFESSIONAL SUCCESS PLANS				
700-020	Success Plan - Professional (Single License) - 1 Year	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, and AI Model Pass for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-021 - Professional - Additional License]</p>	Subscription	\$1,850
700-021	Success Plan - Professional (Additional License) - 1 Year	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, and New AI Model Pass for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-020 - Professional (Single License)]</p>	Subscription	\$1,095
700-026	Success Plan - Professional (Single License) - 3 Years	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, and AI Model Pass for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-027 - Professional - Additional License]</p>	Subscription	\$4,440
700-027	Success Plan - Professional (Additional License) - 3 Years	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, and New AI Model Pass for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-026 - Professional (Single License)]</p>	Subscription	\$2,628

Image-Pro v11.1

ULTIMATE SUCCESS PLANS				
700-022	Success Plan - Ultimate (Single License) - 1 Year	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-023 - Ultimate - Additional License]</p>	Subscription	\$3,895
700-023	Success Plan - Ultimate (Additional License) - 1 Year	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, New AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Supported" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-022 - Ultimate (Single License)]</p>	Subscription	\$1,250
700-028	Success Plan - Ultimate (Single License) - 3 Years	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-029 - Ultimate - Additional License]</p>	Subscription	\$9,348
700-029	Success Plan - Ultimate (Additional License) - 3 Years	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, New AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Supported" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-028 - Ultimate (Single License)]</p>	Subscription	\$3,000

Image-Pro v11.1

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
--------------	--------------	-------------	-------------------------------	------------------

All Products have the option of being delivered with a USB License Dongle, listed as **Physical (P)** or delivered Electronically via **E-Delivery (E)** to be linked to the Computer for the **same price**.

UPGRADES

IMAGE PRO - UPGRADE OPTIONS

100-214	Image-Pro v11.1 - Upgrade from Basic & Insight Products	Upgrade to the latest version of Image-Pro from any previous version of a "Basic" or "Insight" branded product to receive a product with similar features and configuration. If the previous configuration included a Capture Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	P -or- E	\$400
100-215	Image-Pro v11.1 - Upgrade from Plus & Analyzer v7 Products	Upgrade to the latest version of Image-Pro from Image-Pro "Plus" or "Analyzer" v7 to receive a product with similar features and configuration. If the previous configuration included a Network License, a Capture Module, or a 3D Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	P -or- E	\$4,200
100-216	Image-Pro v11.1 - Upgrade from v9, 10, and AutoQuant Products	Upgrade to the latest version of Image-Pro from Image-Pro "Premier" v9 or Image-Pro v10 or even any version of AutoQuant to receive a product with similar features and configuration. If the previous configuration included a Network License, a Capture Module, or a 3D Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	P -or- E	\$2,350
100-217	Image-Pro v11.1 - Upgrade from v11.0	Upgrade to the latest version of Image-Pro from Image-Pro v11.0 to receive a product with similar features, modules, and configuration. If the previous configuration included a License Server, the upgrade to the latest version will also include the new equivalent at no additional charge.	P -or- E	\$1,250
100-218	Image-Pro v11.1 - Upgrade from v9 or v10 Products - with 3D module	Upgrade to the latest version of Image-Pro from Image-Pro "Premier 3D" v9 or Image-Pro v10 with the "3D module" to receive a product with similar features and configuration. If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	P -or- E	\$3,350

Image-Pro v11.1

Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
--------------	--------------	-------------	-------------------------------	------------------

All Products have the option of being delivered with a USB License Dongle, listed as **Physical (P)** or delivered Electronically via **E-Delivery (E)** to be linked to the Computer for the **same price**.

COMPLETE PRICE LIST

PACKAGES

PACKAGES - GENERAL				
Product Code	Product Name	Description	Physical (P) / E-Delivery (E)	List Price (USD)
100-204	Image-Pro v11.1 - Capture Module	Capture images & videos from 150+ devices and then process, adjust, and view images with this powerful and streamlined capture package. INCLUDES: Base Product, 2D Capture Module	P -or- E	\$945
100-205	Image-Pro v11.1 - Basic Package	Extract, compile, review, and report on precise manual measurement data with entry-level Basic bundle. INCLUDES: Base Product, 2D Measurements Module	P -or- E	\$1,575
100-206	Image-Pro v11.1 - 2D Analysis Package	Machine Learning powered segmentation and classification enables precise analysis and automated measurement with the 2D Analysis package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module	P -or- E	\$5,795
100-219	Image-Pro v11.1 - 3D Analysis Package	Visualize stunning 3D views with robust 3D controls and perform precise manual 3D measurements and automated 3D segmentation with the 3D Analysis package. INCLUDES: Base Product, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$7,065
100-207	Image-Pro v11.1 - 2D/3D Analysis Package	Machine Learning powered segmentation and classification enables precise analysis and automated measurement from 2D or 3D images with the 2D/3D Analysis package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$12,545
100-220	Image-Pro v11.1 - AI Complete Analysis Package (2D/3D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D Life Science and Materials applications and is combined with 3D Visualization and Analysis in the AI Complete Analysis Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Protocols - Cell Biology, Protocols - Cell Biology Plus, Protocols - Cell Biology EM, Protocols - Materials, Proocols - Materials EM, AI Deep Learning Extension, Pre-Trained Models - Life Science, Pre-Trained Models - Materials, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$19,520

Image-Pro v11.1

PACKAGES - LIFE SCIENCE				
100-208	Image-Pro v11.1 - Multi-Channel Capture Package	Expand any microscope into a multi-channel fluorescence capture system with this bundle. Control lightsource devices, capture mixed sets of fluorescence and brightfield, and deconvolve using AutoQuant Real-Time Deconvolution, resulting in a crisp, high fidelity fluorescent image set. INCLUDES: Base Product, 2D Capture - Module, Multi-Channel Capture - Extension, & Real-Time Deconvolution - Extension	P -or- E	\$1,665
100-209	Image-Pro v11.1 - Deconvolution Station Package	Import, deconvolve, visualize, and generate animations in a standalone bundle designed to process and visualize fluorescent 3D image stacks. INCLUDES: Base Product, 3D Visualization Module, & AutoQuant Deconvolution Module	P -or- E	\$4,905
100-210	Image-Pro v11.1 - Life Science Package (2D)	Machine Learning powered segmentation and classification enables precise analysis and automated measurement from 2D images with the ML Life Science Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Cell Biology, Cell Biology Plus Protocol Extensions	P -or- E	\$6,875
100-211	Image-Pro v11.1 - AI Life Science Package (2D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Life Science applications with the AI Life Science Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Protocols - Cell Biology, Protocols - Cell Biology Plus, Protocols - Cell Biology EM, AI Deep Learning Extension, Pre-Trained Models - Life Science	P -or- E	\$11,060
100-221	Image-Pro v11.1 - AI Life Science Package (2D/3D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Life Science applications and is combined with 3D Visualization and Analysis in the AI Life Science Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Protocols - Cell Biology, Protocols - Cell Biology Plus, Protocols - Cell Biology EM, AI Deep Learning Extension, Pre-Trained Models - Life Science, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$17,810

Image-Pro v11.1

PACKAGES - MATERIALS				
100-212	Image-Pro v11.1 - Materials Package (2D)	Machine Learning powered segmentation and classification enables precise analysis and automated measurement from 2D images with the ML Materials Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlaysis Module, Materials Protocol Extension	P -or- E	\$6,425
100-213	Image-Pro v11.1 - AI Materials Package (2D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Materials applications with the AI Materials Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlaysis Module, Materials Protocols, & Materials EM Protocol Extension, AI Deep Learning Module, Pre-Trained Materials Models Extension	P -or- E	\$10,655
100-222	Image-Pro v11.1 - AI Materials Package (2D/3D)	Next generation AI Deep Learning technology powers remarkably accurate image segmentation made available as Pre-Trained Models to be combined with efficient analysis protocols for automated analysis in 2D for Materials applications and is combined with 3D Visualization and Analysis in the AI Materials Package. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlaysis Module, Materials Protocols, & Materials EM Protocols Extension, AI Deep Learning Module, Pre-Trained Materials Models Extension, 3D Visualization Module, & 3D Measure & Analyze Extension	P -or- E	\$17,405

Image-Pro v11.1

INDIVIDUAL

BASE PRODUCT				
100-201	Image-Pro - v11.1 Base	The Image-Pro Base is required for all modules and extensions. Includes File Management, 2D Viewer, Process/Adjust, Usability, Data Management, Reports, Auditing, Scripting LINK TARGETS: [Windows PC, Green USB dongle (optional), or License Server Module (sold separately); Purple USB dongles not supported]	P -or- E	\$350

MODULES & EXTENSIONS				
300-024	2D Measurements - v11.1 Module	Extract precise 2D manual measurement data using semi-automated tools for fast and accurate results. Includes Direct & Assisted Measurements, Auto Shape & Edge Detection, Classification, & Color Measurement. REQUIRES: [100-201 - Image-Pro - Base]	P -or- E	\$1,400
300-025	2D Automated Analysis - v11.1 Module	Tackle complex analysis across multiple images/channels with Protocols or Free-Form analysis with ML-based and optional DL-based tools. Includes Smart Segmentation, Classification, 2D Tracking, & the Essentials Protocol Collection with 3 Protocols. REQUIRES: [100-201 - Image-Pro - Base]	P -or- E	\$4,690
300-026	Protocols - Cell Biology - v11.1 Extension	Analyze complex Cell Biology-based measurements with an efficient and repeatable process. Can be enhanced with AI Deep Learning Extension. Includes: Cell Apoptosis, Cell Autophagy, Cell Confluence, Cell Count, Cell Morphology, Cell Proliferation, Cell Transfection, Lipid Droplet, Live/Dead Cell Viability, & Ring Region. REQUIRES: [300-025 - 2D Automated Analysis - Module]	P -or- E	\$550
300-027	Protocols - Cell Biology Plus - v11.1 Extension	Analyze advanced Cell Biology-based measurements with an efficient and repeatable process. Includes: Angiogenesis, Colocalization, Neurite Outgrowth, Translocation, & Wound Healing. REQUIRES: [300-025 - 2D Automated Analysis - Module]	P -or- E	\$650
300-028	Protocols - Materials - v11.1 Extension	Analyze complex Materials-based measurements with an efficient and repeatable process. Includes: Composite Materials, Fiber Separation, Fiber Thickness, Grain, Layers, Particles Phase, Particle Size, Phases, & Pores. REQUIRES: [300-025 - 2D Automated Analysis - Module]	P -or- E	\$700

Image-Pro v11.1

300-032	Protocols - Cell Biology EM - v11.1 Extension	Analyze complex EM Cell Biology-based measurements with an efficient and repeatable process. Includes: Mitochondria (EM), Nuclei and Chromatin (EM), Nuclei and Mitochondria (EM), & Nuclei (EM). REQUIRES: [300-025 - 2D Automated Analysis - Module]	P -or- E	\$650
300-033	Protocols - Materials EM - v11.1 Extension	Analyze complex EM Materials-based measurements with an efficient and repeatable process. Includes: Composite Materials (EM), Fiber Separation (EM), Fiber Thickness (EM), Grain (EM), Layers (EM), Particle Phase (EM), Particle Size (EM), Phases (EM), & Pores (EM). REQUIRES: [300-025 - 2D Automated Analysis - Module]	P -or- E	\$700
300-029	AI Deep Learning - v11.1 Extension	Next generation AI Deep Learning technology for remarkably accurate image segmentation. Use stand-alone or within Analysis Protocols. Use Pre-Trained Models (Sold Separately) or train your own models with the intuitive Model Trainer. Includes "General Objects" Pre-Trained Model. REQUIRES: [300-025 - 2D Automated Analysis - Module]	P -or- E	\$3,500
300-030	Pre-Trained Models - Life Science - v11.1 Extension	Pre-Trained Models for Life Science-based samples. Use to instantly detect and measure or as a starting point for fine-tuning a duplicate model. Includes: Fluorescent Cells, Dendrites, Nuclei, & Overlapping Nuclei, Label Free Cells, & Nuclei, Hematoxylin Nuclei, TEM Mitochondria, TEM Nuclei, and Retinal Blood Vessels. REQUIRES: [300-029 - 2D Automated Analysis - Module]	P -or- E	\$500
300-031	Pre-Trained Models - Materials - v11.1 Extension	Pre-Trained Models for Materials-based samples. Use to instantly detect and measure or as a starting point for fine-tuning a duplicate model. Includes: Versatile Grains, Versatile Particles, Ti6Al4V Titanium Alloy, Ti6Al4V Alpha Phase, and Widmanstätten Patterns. REQUIRES: [300-029 - 2D Automated Analysis - Module]	P -or- E	\$500
300-034	2D Capture - v11.1 Module	Capture images & videos from 150+ devices, perform live measurement, and live processing with ease. Includes the Live Image Workspace, Instant & HDR Capture, Recording of Time-Lapse Events, and basic Presets, Enhancements, and Preview/Acquire settings. REQUIRES: [100-201 - Image-Pro - Base]	P -or- E	\$700
300-035	Real-Time Deconvolution - v11.1 Extension	Elevate any capture system with Real-Time Deconvolution to deliver crisp, high fidelity, fluorescent images during preview and capture. Includes 2D iterative Deconvolution & Wiener Filter algorithms. REQUIRES: [300-034 - 2D Capture - Module]	P -or- E	\$400
300-036	Multi Channel Capture - v11.1 Extension	Expand any microscope with Multi-Channel Capture making it capable of capturing multiple channels into an image set. Includes LightSource Device Control, Acquisition of mixed sets with Fluorescence and Brightfield, and Auto & Manual capture modes. REQUIRES: [300-034 - 2D Capture - Module]	P -or- E	\$400

Image-Pro v11.1

300-037	3D Visualization - v11.1 Module	<p>Visualize 3D image stacks with full control over appearance for the generating of stunning animations with complex movements & renderings. Includes 3D viewer with advanced controls, Min/Max projections, clipping planes, slicer viewing, & Movie Maker.</p> <p>REQUIRES: [100-201 - Image-Pro - Base]</p>	P -or- E	\$2,000
300-038	3D Measure & Analyze - v11.1 Extension	<p>Measure & Analyze with familiar tools across multiple dimensions. Includes iso-surface based object segmentation, 3d volume measurements, guided segmentation, 3D colocalization, and manual point-to-point distances, angles, and more in 3D space.</p> <p>REQUIRES: [300-037 - 3D Visualization - Module]</p>	P -or- E	\$5,500
300-039	3D Big Data Visualization - v11.1 Extension	<p>Convert large datasets captured by high-volume imaging techniques to multi-resolution pyramidal formats for quick and easy visualization. Includes background file conversion, terabyte-sized file support, and multi-resolution file visualization.</p> <p>REQUIRES: [300-037 - 3D Visualization - Module]</p>	P -or- E	\$2,000
300-040	AutoQuant Deconvolution - v11.1 Module	<p>Restore image fidelity, enhance quality, and reverse distortions introduced during the image acquisition process. Includes 2D & 3D Algorithms, Theoretical & Measured PSF models, Spherical Abberation correction, GPU acceleration, and Batch processing.</p> <p>REQUIRES: [100-201 - Image-Pro - Base]</p>	P -or- E	\$3,100
300-023	License Server - v11.1 Module	<p>Distribute software licenses to teams with a stand-alone License Server and remove USB dongle and Node-Locked licenses from the network. Includes floating or reserved to workstation options, flexible bundle assignment, auto-license return after user time outs, options for VPN and Virtual Machine support.</p> <p>Includes Live License Server Onboarding</p> <p>LINK TARGETS: [Windows, Windows Server, Linux, VM's ; See IT Checklist for complete list]</p>	E	\$2,500

Image-Pro v11.1

UPGRADES

IMAGE PRO - UPGRADE OPTIONS				
100-214	Image-Pro v11.1 - Upgrade from Basic & Insight Products	<p>Upgrade to the latest version of Image-Pro from any previous version of a "Basic" or "Insight" branded product to receive a product with similar features and configuration.</p> <p>If the previous configuration included a Capture Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.</p>	P -or- E	\$400
100-215	Image-Pro v11.1 - Upgrade from Plus & Analyzer v7 Products	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Plus" or "Analyzer" v7 to receive a product with similar features and configuration.</p> <p>If the previous configuration included a Network License, a Capture Module, or a 3D Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.</p>	P -or- E	\$4,200
100-216	Image-Pro v11.1 - Upgrade from v9, 10, and AutoQuant Products	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Premier" v9 or Image-Pro v10 or even any version of AutoQuant to receive a product with similar features and configuration.</p> <p>If the previous configuration included a Network License, a Capture Module, or a 3D Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.</p>	P -or- E	\$2,350
100-217	Image-Pro v11.1 - Upgrade from v11.0	<p>Upgrade to the latest version of Image-Pro from Image-Pro v11.0 to receive a product with similar features, modules, and configuration.</p> <p>If the previous configuration included a License Server, the upgrade to the latest version will also include the new equivalent at no additional charge.</p>	P -or- E	\$1,250
100-218	Image-Pro v11.1 - Upgrade from v9 or v10 Products - with 3D module	<p>Upgrade to the latest version of Image-Pro from Image-Pro "Premier 3D" v9 or Image-Pro v10 with the "3D module" to receive a product with similar features and configuration.</p> <p>If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.</p>	P -or- E	\$3,350

ACCESSORIES

ACCESSORIES				
700-016	USB Dongle	A single secure USB dongle (green) with license linking capability for Image-Pro. Use the License Manager to select this USB dongle as a target for linking an Image-Pro v11 or higher license.	P	\$700
700-017	USB Dongle - Demo NFR - Black	A single secure USB dongle (black) that can only be linked to a Not For Resale (NFR) demo license. Use this dongle for Demo Purposes ONLY. Use the License Manager to select this USB dongle as a target for linking an Image-Pro v11 or higher NFR Demo License.	P	\$150

Image-Pro v11.1

Image-Pro v11.1

SERVICE - PLANS

STANDARD SUCCESS PLANS				
700-018	Success Plan - Standard (Single License) - 1 Year	<p>This plan provides a single onboarding session and continuous software support for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-019 - Standard - Additional License]</p>	Subscription	\$1,450
700-019	Success Plan - Standard (Additional License) - 1 Year	<p>Add an additional license to your Plan and receive continuous software support for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>REQUIRES: [700-018 - Standard (Single License)]</p>	Subscription	\$150
700-024	Success Plan - Standard (Single License) - 3 Years	<p>This plan provides a single onboarding session and continuous software support for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-025 - Standard - Additional License]</p>	Subscription	\$3,480
700-025	Success Plan - Standard (Additional License) - 3 Years	<p>Add an additional license to your Plan and receive continuous software support for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue suport services.</p> <p>REQUIRES: [700-024 - Standard (Single License)]</p>	Subscription	\$360

Image-Pro v11.1

PROFESSIONAL SUCCESS PLANS				
700-020	Success Plan - Professional (Single License) - 1 Year	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, and AI Model Pass for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-021 - Professional - Additional License]</p>	Subscription	\$1,850
700-021	Success Plan - Professional (Additional License) - 1 Year	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, and New AI Model Pass for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-020 - Professional (Single License)]</p>	Subscription	\$1,095
700-026	Success Plan - Professional (Single License) - 3 Years	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, and AI Model Pass for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-027 - Professional - Additional License]</p>	Subscription	\$4,440
700-027	Success Plan - Professional (Additional License) - 3 Years	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, and New AI Model Pass for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-026 - Professional (Single License)]</p>	Subscription	\$2,628

Image-Pro v11.1

ULTIMATE SUCCESS PLANS				
700-022	Success Plan - Ultimate (Single License) - 1 Year	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-023 - Ultimate - Additional License]</p>	Subscription	\$3,895
700-023	Success Plan - Ultimate (Additional License) - 1 Year	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, New AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 1 year from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Supported" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-022 - Ultimate (Single License)]</p>	Subscription	\$1,250
700-028	Success Plan - Ultimate (Single License) - 3 Years	<p>This plan provides a single onboarding session, continuous software support, Software Upgrade Protection, the Protocol Pass, AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Active" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>FOR ADDITIONAL LICENSES USE: [700-029 - Ultimate - Additional License]</p>	Subscription	\$9,348
700-029	Success Plan - Ultimate (Additional License) - 3 Years	<p>Add an additional license to your Plan and receive continuous software support, Software Upgrade Protection, the New Protocol Pass, New AI Model Pass, Uptime Guarantee, License Server Support, Training, and Application Services for a term of 3 years from the time of purchase. Services are provided by Media Cybernetics professional support staff by email, phone, or screenshare to help keep your software running at peak performance year-over-year.</p> <p>COVERAGE: A single software license listed as "Supported" by the Software Life Cycle Policy. RENEWAL: Annual renewal required to continue support services.</p> <p>REQUIRES: [700-028 - Ultimate (Single License)]</p>	Subscription	\$3,000