



Image-Pro v11 - Summary - Page 1

Product Code	Product Name	Requires an existing...	Includes	Delivery	Pricing (USD)
Base Product					
100-001	Image-Pro - Base			Electronic	\$300
Modules & Extensions					
300-011	2D Measurements - Module	Base (100-001)		Electronic	\$1,400
300-012	2D Automated Analysis - Module	Base (100-001)		Electronic	\$4,690
300-013	Protocols - Cell Biology - Extension	2D Automated Analysis Module (300-012)		Electronic	\$500
300-014	Protocols - Cell Biology Plus - Extension	2D Automated Analysis Module (300-012)		Electronic	\$600
300-015	Protocols - Materials - Extension	2D Automated Analysis Module (300-012)		Electronic	\$650
300-016	2D Capture - Module	Base (100-001)		Electronic	\$700
300-017	Real-Time Deconvolution - Extension	2D Capture Module (300-016)		Electronic	\$400
300-018	Multi-Channel Capture - Extension	2D Capture Module (300-016)		Electronic	\$400
300-019	3D Visualization - Module	Base (100-001)		Electronic	\$2,000
300-020	3D Measure & Analyze - Extension	3D Visualization Module (300-019)		Electronic	\$5,500
300-021	3D Big Data Visualization - Extension	3D Visualization Module (300-019)		Electronic	\$2,000
300-022	AutoQuant Deconvolution - Module	Base (100-001)		Electronic	\$3,000
300-004	License Server - Module	Base (100-001)		Electronic	\$2,000
Bundles - General					
100-101	Image-Pro - Basic Bundle		100-001, 300-011	Electronic	\$1,530
100-102	Image-Pro - 2D Analysis Bundle		100-001, 300-011, 300-012	Electronic	\$5,750
100-103	Image-Pro - 2D/3D Analysis Bundle		100-001, 300-011, 300-012, 300-019, 300-020	Electronic	\$12,500
Bundles - Life Science					
100-104	Image-Pro - Deconvolution Station Bundle		100-001, 300-019, 300-022	Electronic	\$4,770
100-105	Image-Pro - 2D Cell Analysis Bundle		100-001, 300-011, 300-012, 300-013, 300-014, 300-022	Electronic	\$9,440
Bundles - Industrial					
100-106	Image-Pro - 2D Materials Analysis Bundle		100-001, 300-011, 300-012, 300-015	Electronic	\$6,336
Upgrades					
100-016	Image-Pro - Upgrade from Basic & Insight Products	Valid Serial Number	100-101, Optional: 300-016	Electronic	\$400
100-017	Image-Pro - Upgrade from Plus & Analyzer v7 Products	Valid Serial Number	100-102 -or- 100-103, Optional: 300-016 & 300-004	Electronic	\$4,200
100-018	Image-Pro - Upgrade from v9, 10, and AutoQuant Products	Valid Serial Number	100-104 -or- 100-102, Optional: 300-016 & 300-004	Electronic	\$2,350
100-019	Image-Pro - Upgrade from v9 or v10 Products - with 3D module	Valid Serial Number	100-103, Optional: 300-021 & 300-016 & 300-004	Electronic	\$3,350
Support & Accessories					
700-002	Upgrade Protection Plan	Base (100-001)		Electronic	\$945
700-016	USB Dongle		Single USB dongle (green) w/ license linking capability	Physical	\$700
700-017	USB Dongle - Demo NFR - Black		Single USB dongle (black) w/ NFR license linking capability	Physical	\$150



Image-Pro v11 - Detailed Descriptions - Pages 2 - 5

Code	Icon	Product Name	Description	Delivery	Pricing (USD)
Base Product					
100-001		Image-Pro - Base	<p>The Image-Pro Base is required for all modules and extensions. Includes File Management, 2D Viewer, Process/Adjust, Usability, Data Management, Reports, Auditing, Scripting</p> <p>LINK TARGETS: [Windows PC, Green USB dongle (sold separately), or License Server Module (sold separately); Purple USB dongles not supported]</p>	Electronic	\$300
Modules & Extensions					
300-011		2D Measurements - Module	<p>Extract precise 2D manual measurement data using semi-automated tools for fast and accurate results. Includes Direct & Assisted Measurements, Auto Shape & Edge Detection, Classification, & Color Measurement.</p> <p>REQUIRES: [100-001 - Image-Pro - Base]</p>	Electronic	\$1,400
300-012		2D Automated Analysis - Module	<p>Tackle complex analysis across multiple images/channels with Protocols or Free-Form analysis with ML-based tools. Includes Smart Segmentation, Classification, 2D Tracking, & the Essentials Protocol Collection with 3 Protocols.</p> <p>REQUIRES: [100-001 - Image-Pro - Base]</p>	Electronic	\$4,690
300-013		Protocols - Cell Biology - Extension	<p>Analyze complex Cell Biology-based measurements with an efficient and repeatable process. Includes 10 Protocols: Apoptosis, Autophagy, Cell Count, Cell Morphology, Cell Proliferation, Confluence, Lipid Droplet, Live/Dead, Ring, & Transfection.</p> <p>REQUIRES: [300-012 - 2D Automated Analysis - Module]</p>	Electronic	\$500
300-014		Protocols - Cell Biology Plus - Extension	<p>Analyze complex Cell Biology-based measurements with an efficient and repeatable process. Includes 5 Protocols: Angiogenesis, Colocalization, Neurite Outgrowth, Translocation, & Wound Healing.</p> <p>REQUIRES: [300-012 - 2D Automated Analysis - Module]</p>	Electronic	\$600
300-015		Protocols - Materials - Extension	<p>Analyze complex Materials-based measurements with an efficient and repeatable process. Includes 10 Protocols: Composite Materials, Layers Concentric, Layers Curved, Layers Parallel, Fiber Separation, Fiber Thickness, Particles, Particles Phase, Phase, & Pores.</p> <p>REQUIRES: [300-012 - 2D Automated Analysis - Module]</p>	Electronic	\$650
300-016		2D Capture - Module	<p>Capture images & videos from 150+ devices, perform live measurement, and live processing with ease. Includes the Live Image Workspace, Instant & HDR Capture, Recording of Time-Lapse Events, and basic Presets, Enhancements, and Preview/Acquire settings.</p> <p>REQUIRES: [100-001 - Image-Pro - Base]</p>	Electronic	\$700
300-017		Real-Time Deconvolution - Extension	<p>Elevate any capture system with Real-Time Deconvolution to deliver crisp, high fidelity, fluorescent images during preview and capture. Includes 2D iterative Deconvolution & Wiener Filter algorithms.</p> <p>REQUIRES: [300-016 - 2D Capture - Module]</p>	Electronic	\$400



Code	Icon	Product Name	Description	Delivery	Pricing (USD)
300-018		Multi-Channel Capture - Extension	Expand any microscope with Multi-Channel Capture making it capable of capturing multiple channels into an image set. Includes LightSource Device Control, Acquisition of mixed sets with Fluorescence and Brightfield, and Auto & Manual capture modes. REQUIRES: [300-016 - 2D Capture - Module]	Electronic	\$400
300-019		3D Visualization - Module	Visualize 3D image stacks with full control over appearance for the generating of stunning animations with complex movements & renderings. Includes 3D viewer with advanced controls, Min/Max projections, clipping planes, slicer viewing, & Movie Maker. REQUIRES: [100-001 - Image-Pro - Base]	Electronic	\$2,000
300-020		3D Measure & Analyze - Extension	Measure & Analyze with familiar tools across multiple dimensions. Includes iso-surface based object segmentation, 3d volume measurements, guided segmentation, 3D colocalization, and manual point-to-point distances, angles, and more in 3D space. REQUIRES: [300-019 - 3D Visualization - Module]	Electronic	\$5,500
300-021		3D Big Data Visualization - Extension	Convert large datasets captured by high-volume imaging techniques to multi-resolution pyramidal formats for quick and easy visualization. Includes background file conversion, terabyte-sized file support, and multi-resolution file visualization. REQUIRES: [300-019 - 3D Visualization - Module]	Electronic	\$2,000
300-022		AutoQuant Deconvolution - Module	Restore image fidelity, enhance quality, and reverse distortions introduced during the image acquisition process. Includes 2D & 3D Algorithms, Theoretical & Measured PSF models, Spherical Abberation correction, GPU acceleration, and Batch processing. REQUIRES: [100-001 - Image-Pro - Base]	Electronic	\$3,000
300-004		License Server - Module	Distribute software licenses to teams with a stand-alone License Server and remove USB dongle and Node-Locked licenses from the network. Includes floating or reserved to workstation options, flexible bundle assignment, auto-license return after user time outs, options for VPN and Virtual Machine support. LINK TARGETS: [Windows, Windows Server, Linux, VM's ; See IT Checklist for complete list]	Electronic	\$2,000



Code	Icon	Product Name	Description	Delivery	Pricing (USD)
Bundles - General					
100-101		Image-Pro - Basic Bundle	Extract, compile, review, and report on precise manual measurement data with entry-level Basic bundle. INCLUDES: Base Product, 2D Measurements Module	Electronic	\$1,530
100-102		Image-Pro - 2D Analysis Bundle	Extract, segment, classify, and analyze precise manual and automated measurement data with the 2D Analysis bundle. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module	Electronic	\$5,750
100-103		Image-Pro - 2D/3D Analysis Bundle	Extract, segment, classify, and analyze precise manual and automated measurement data from 2D or 3D images with the 2D/3D Analysis bundle. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, 3D Visualization Module, & 3D Measure & Analyze Extension	Electronic	\$12,500
Bundles - Life Science					
100-104		Image-Pro - Deconvolution Station Bundle	Import, deconvolve, visualize, and generate animations in a standalone bundle designed to process and visualize fluorescent 3D image stacks. INCLUDES: Base Product, 3D Visualization Module, & AutoQuant Deconvolution Module	Electronic	\$4,770
100-105		Image-Pro - 2D Cell Analysis Bundle	Extract, segment, classify, analyze, and report precise 2D manual and automated cell measurement data with the 2D Cell Analysis bundle. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Analysis Module, Cell Biology & Cell Biology Plus Protocol Extensions, & AutoQuant Deconvolution Module	Electronic	\$9,440
Bundles - Industrial					
100-106		Image-Pro - 2D Materials Analysis Bundle	Extract, segment, classify, analyze, and report precise 2D manual and automated materials measurement data with the 2D Materials Analysis bundle. INCLUDES: Base Product, 2D Measurements Module, 2D Automated Anlaysis Module, Materials Protocol Extension	Electronic	\$6,336








Code	Icon	Product Name	Description	Delivery	Pricing (USD)
Upgrades					
100-016		Image-Pro - Upgrade from Basic & Insight Products	Upgrade to the latest version of Image-Pro from any previous version of a "Basic" or "Insight" branded product to receive a product with similar features and configuration. If the previous configuration included a Capture Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	Electronic	\$400
100-017		Image-Pro - Upgrade from Plus & Analyzer v7 Products	Upgrade to the latest version of Image-Pro from Image-Pro "Plus" or "Analyzer" v7 to receive a product with similar features and configuration. If the previous configuration included a Network License, a Capture Module, or a 3D Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	Electronic	\$4,200
100-018		Image-Pro - Upgrade from v9, 10, and AutoQuant Products	Upgrade to the latest version of Image-Pro from Image-Pro "Premier" v9 or Image-Pro v10 or even any version of AutoQuant to receive a product with similar features and configuration. If the previous configuration included a Network License, a Capture Module, or a 3D Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	Electronic	\$2,310
100-019		Image-Pro - Upgrade from v9 or v10 Products - with 3D module	Upgrade to the latest version of Image-Pro from Image-Pro "Premier 3D" v9 or Image-Pro v10 with the "3D module" to receive a product with similar features and configuration. If the previous configuration included a Network License or a Capture Module, the upgrade to the latest version will also include the new equivalent modules at no additional charge.	Electronic	\$3,350
Support & Accessories					
700-002		Upgrade Protection Plan	Upgrade Protection ensures your product is always up-to-date with the latest features. While under coverage, all upgrades to features within your product configuration will be provided at no additional cost. Annual Renewal required.	Electronic	\$945
700-016		USB Dongle	A single secure USB dongle (green) with license linking capability for Image-Pro. Use the License Manager to select this USB dongle as a target for linking an Image-Pro v11 or higher license.	Physical	\$700
700-017		USB Dongle - Demo NFR - Black	A single secure USB dongle (black) that can only be linked to a Not For Resale (NFR) demo license. Use this dongle for Demo Purposes ONLY. Use the License Manager to select this USB dongle as a target for linking an Image-Pro v11 or higher NFR Demo License.	Physical	\$150



Image-Pro Plus & Analyzer - Page 6

Product Code	Product Description	Includes	Delivery	Pricing (USD)
--------------	---------------------	----------	----------	---------------

NEW

Single Permanent Licenses				
102-001	Image Pro Plus v7 - Single User License - USB Dongle	Installion USB Flash Drive & USB License Dongle	Physical	\$7,650
102-005	Image-Pro Analyzer v7 - Single User License - USB Dongle	Installion USB Flash Drive & USB License Dongle	Physical	\$6,480

Accessories				
700-015	Replacement USB Dongle for Image-Pro v7 - v10, & AutoQuant X3	A single secure USB dongle	Physical	\$675

Retired Part Numbers

As of February 28th, 2023



Retired Part Numbers - Pages 7 - 10

Product Code	Product Description	Pricing (USD)	Reason
IMAGE-PRO PLUS/ANALYZER NEW			
Single Permanent Licenses			
102-002	Image-Pro Plus DVD Option	\$330	No longer offering
102-003	Image-Pro Plus User Guide	\$550	No longer offering
Upgrades			
102-004	Image-Pro Plus Upgrade to Version 7	\$4,215	Only new sales of v7 or upgrades to v11
102-006	Image-Pro Analyzer Upgrade to Version 7	\$3,635	No longer offering
Maintenance and Support			
700-014	Lost, Broken, & Replacement Key for Image-Pro Family Version 7 and older	\$675	Replaced by 700-015



Retired Part Numbers - Pages 7 - 10

IMAGE-PRO v10 NEW

Single Permanent Licenses		
100-002	Image-Pro Single User License (USB)	\$5,775
100-006	Image-Pro DVD option (requires 100-002)	\$315
901-003	Image-Pro Basic Configuration - Permanent Single License - USB Dongle	\$1,570
901-004	Image-Pro Basic - USB - 10 Pack Pricing, Must Order Qty of 10	\$1,495
Multi-Seat Network Server License		
100-009	Network Server License (Soft License) includes One Image-Pro Single User License Seat	\$7,350
100-010	Image-Pro Single User License (Software License) for Network Server License (Requires 100-009 or 100-009 Serial Number)	\$5,775
Modules		
300-001	Capture Module for Image-Pro Platform	\$1,260
300-003	3D Module for Image-Pro Platform	\$7,875
300-005	3D Module for Image-Pro Platform with Big Data Visualization	\$9,975
Apps		
601-001	Inter-Metallic Corrosion (IMC) APP	\$5,250
601-002	Grains APP	\$1,575
601-003	Fiber Thickness APP	\$2,625

v10 replaced by v11
No longer offering
Replaced by 100-001 Basic Bundle
Replaced by 100-001 Basic Bundle

v10 replaced by v11 License Server
v10 replaced by v11 License Server

v10 replaced by v11
v10 replaced by v11
v10 replaced by v11

Now a Free App
Now a Free App
Now a Free App & Paid Protocol in v11

IMAGE-PRO v10
UPGRADES

Single Permanent License Upgrades		
100-003	Image-Pro Single User License Upgrade (Image-Pro Premier, Image-Pro Premier 3D, Image-Pro Insight)	\$2,310
100-007	Image-Pro Basic Single User License Upgrade	\$4,200
100-008	Image-Pro Legacy Single User License Upgrade (Includes Dongle) Versions (Image-Pro Plus and Image-Pro Analyzer)	\$3,885
300-010	Image-Pro 3D Module Big Data Visualization Upgrade (Requires 300-003)	\$3,150
Network Server License Upgrades		
100-011	Image-Pro Network License Server Incudes Image-Pro Single User License Seat	\$3,675
100-012	Image-Pro Network License Server Seat Upgrade (Requires 100-011)	\$3,150

v10 replaced by v11
v10 replaced by v11
v10 replaced by v11
v10 replaced by v11

v10 replaced by v11 License Server
v10 replaced by v11 License Server

IMAGE-PRO v10 SUPPORT LICENSES

Maintenance and Support		
700-001	Lost, Broken or Replacement Key	\$675

Replaced by 700-015

IMAGE-PRO v10 NFR PARTNER DEMO LICENSES

Not For Resale (NFR) Licenses		
100-004	Image-Pro Demo License (NFR)	\$110
100-005	Image-Pro NFR Key Renewal	\$0

Replaced by 700-017
No longer valid



Retired Part Numbers - Pages 7 - 10

AUTOQUANT NEW

Single Permanent Licenses		
AQX3-CWF-USB	AutoQuant X3 - Confocal & Widefield	\$7,000
AQX3-CWF-ADD-USB	AutoQuant X3 - Confocal & Widefield - Permanent Single License - 1 Additional Copy Discount - USB Dongle	\$5,600
AQX-GPU-CWF	2D & 3D Confocal and Widefield Blind and Non-Blind deconvolution, AutoQuant X3 Confocal & Widefield GPU Processing Module. Requires the latest version of AutoQuant X3-Confocal & Widefield version 3.1	\$2,500

Replaced by v11 AutoQuant Module

Replaced by v11 AutoQuant Module

Replaced by v11 AutoQuant Module

AUTOQUANT UPGRADES

Single Permanent Licenses Upgrades		
AQX3-CWF-U1	Upgrade to AQX3 CWF from X2 CWF	\$1,750
AQX3-CWF-U2	Upgrade to AQX3 CWF from X CWF and older	\$2,450

Upgrade to v11 using 100-018

Upgrade to v11 using 100-018

AUTOQUANT SUPPORT LICENSES

Replacement Keys - Single Licenses		
AQX-RKEY	AutoQuant Family Exchange USB Key, Under Warranty	\$0
AQX-BKEY	AutoQuant Family Exchange USB Key, Expired Warranty	\$250
AQX-LKEY	AutoQuant Family Replacement USB Key	\$500

Replaced by 700-017

Replaced by 700-017

Replaced by 700-017

AUTOQUANT NFR PARTNER DEMO LICENSES

Not For Resale (NFR) Licenses		
AQX3-CWF-NFR-USB	AutoQuant X3 – Total - Not For Resale	\$150

Replaced by v11 using 700-017